# Sling assessment and how to meet the requirements in the ISO standard

# SAFER Handling Conference BRISBANE 27th of May 2014





### Ole Lindahl Thy

Education Manager, Etac AB

I live in Lillehammer, Norway. (*Not* Sweden) This is my wife in our "backyard"...



Education



### Introduction

Choosing the right sling is essential to secure a safe and comfortable transfer for both the caregiver and the client.

It is important to make an individual assessment in every situation where a client is mechanical lifted.



Education



## Introduction, - Causes of hoisting incidents

Cause of incident	Number	%	
Fall from hoist	99	61	
Equipment failure	31	19	
Incorrect equipment used	8	5	
Left unattended	7	4	
Hoist overturn	5	3	
Others	13	8	
Total	163	100	

Reported incidents (2001-2007) to the Health & safety executive (HSE) in the UK, in which people have been injured while being moved using hoisting equipment **(Boulton 2009)** 



## **Causes of incidents**

Fall from hoist

- Use of a sling that is too large, no adequate support, a large aperture
- Using toileting sling for all hoisting task, leading to lack of support
- Failing to attach or attach securely one or more loops/clips

### Equipment failure

- Slings with damaged stitching, loops and clips
- "worn out" hoist

Incorrect equipment used

 Wrong size of sling, caused by appropriate sling was not being available

Education

Use of standing hoist, unsuitable for the patient



# The requirements in the ISO Standard

The international standard for hoists and slings is the ISO 10535-2006. Most manufacturers of hoists and slings are testing their products according to this standard.





# The requirements for body-support units

ISO 10535-2006

- 4 General requirements and test methods
- 4.4 Requirements for body-support units

The manufacturer of the body-support unit shall indicate with which hoist(s) and spreader bar(s) it is compatible in order to ensure a safe combination





## A combinationlist for spreaderbar & slings

The manufacturer of the body-support unit shall indicate with which hoist(s) and spreader bar(s) it is compatible in order to ensure a safe combination

Combination list				Suspension Rgo Sling Medium back						
Suspension		Rgo Sl	ing Mediu	m back		Smart		XXS - XS	S - L	XL - XXI
Mover 180		/ Rgo	Sling Hig	h back		9	4-pt			
Contract 100	2-pt Small	<ul> <li>Image: A state of the state of</li></ul>	5-1	AL - AAL	Ø.		Medium 440mm	~	~	~
<b>1</b>	2-pt Medium 440mm	~	~			Nomad				
a constant	2-pt Large 540mm		~		9		4-pt Medium	~	~	~
a ser	4-pt Medium 440mm	~	~	~		Nova	2-pt			
	4-pt Large 540mm		~	~	6		Small 2-pt Medium	<ul><li>✓</li></ul>	~	
Mover 205	-				(	1	2-pt		~	
	2-pt Small 340mm	~			Œ		4-pt Medium		~	~
	2-pt Medium 440mm	~	~		Æ		4-pt		-	~
	2-pt Large 540mm		~		Slinga	U uide	Laige			
	4-pt Medium 440mm	~	~	~	The size weight	e to choose d , function and afe Working I	epends pa body circ oad) 300	rtly on t umferer ka	he use ice/size	er's e.
	4-pt Large 540mm		~	~	Size XXL	Weight (Kg) 230-300				
Mover 300					XL	160-240	1			
	4-pt Medium 440mm	~	~	~	L M	90-160 45-95				
	4-pt Large 540mm		~	~	s XS XXS	25-50 17-25 12-17				
Partner 255										
0	4-pt Medium 440mm	~	~	~						
	4-pt Large 540mm		~	~						
-				<u> </u>			Educ	atio	n (	×



## Non-compatible spreaderbars & slings

The connection points of the spreader-bar and the sling are compatible

Example of non-compatible connection points:





## Requirements for risk assessment

ISO 10535-2006

- 8 Non-rigid body-support units specific requirements and test methods
- 8.4.2 Instruction for use

A statement shall be included, warning the user that a risk assessment shall be carried out to ensure that the correct size, type and shape of body-support unit is being used for the patient



### Information required on the sling or in the User Manual



### Information on the sling or in the User Manual

### The demands in ISO 10535 – 2006

Information that are allowed given in the instruction of use:

- Field of application, direction of use
- The method of lifting
- Materials used
- A warning not to use a damaged or badly worn sling
- The method by which the sling can be removed
- A statement telling the carer that a risk assessment shall be carried out for each patient.







Never use a faulty or damaged sling as it can break and cause personal injury. Destroy and discard damaged and old slings



C14. Remove the log support by folding the leg support under the users thigh and pull out.
C15. Stand new too the chair and lean the user forward while supporting with other hand, and hat the slarg supporting with other hand, and hat the slarg support of the the chair if you pull the slarg to hard, the user may fall for wards and be injured. Always hold one hand around the user.

It is important that the sling has been tested with the individual user and for the intended lifting situation. Make a decision on whether one or more assistants are required.

Plan the lifting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.

Assess the risks and take notes. You as a carer are responsible for the safety of the user.



C16. Release the loops from the suspension. Turn the user towards you until herbin is hying safely on their side. Fold the sing in the middle and position it behind the user's back. Then turn the user over on the opposite side and pull the sing away. Be careful when repositioning the user over on the other side. Make sure herbin does not filp over and rolls out of bad





### Periodic Inspection Iso 10535 – 2006, Annex B: 8 -11

**B.8** Periodic inspection of the non-rigid body-support unit should be undertaken at the time intervals stated by the manufacturer, but at least every 6 months. More frequent inspections may be required where a non-rigid body-support unit is used or cleaned more frequently than normal.

**B.9** Inspections should be performed by a person who is suitably and properly qualified and well acquainted with the design, use and care of the body-support unit.

The inspection should be to find signs of damage, wear or potential failure.

- B.10 The inspection record should be retained safely for examination in the event of an incident.
- B.11 The inspection record should include the following information:
- date of inspection;
- identification details and serial number of the body-support unit;
- information about the condition of the body-support unit;
- date next inspection is due;
- identification and signature of the inspector.





### RgoSling Periodic Inspection

### Scope

Periodic inspection must be carried out at least every 6 months. More regular inspection may be required if the sling is used or washed more often than is normal, inspection shall be performed. by trained personell. Use this checklist to find signs of damage, wear or potential failure. Verify that the sling is safe to use until next periodic inspection. Inspection record should retained safely for examination in case of an accident.

Prouct Information (use list if more then one): Sling Model:

### (Ex. Toilet lowback)

Serial no.:

Art.No: DXXS DXS DS DM Size: 

that does not pass inspection.

Destroy and discard damaged and old slings

### Checklist

OK	Fallure	N/A		
			1. Product label It is possible to read; Maxioad (SWL), Article no., Serial number, Size, Was	hing symbols.
			2. Fabric No tears or frays, burnmarks, holes, color variations, bad smell or other sit	ins of wear
			3. Padding No shrinkaga, bulks, bolds or folds.	<u> </u>
			4. Seams No loose threads, or damaged seams/stiches.	000
			5. Lifting Straps and edge band No tears or frays, stretchmarks, burnmarks or other signs of wear	() ()
			<ol> <li>Guide handles and pockets</li> <li>No tear, frays or locse stitches.</li> </ol>	_0-
			7. Details Ladderlocks, buckles or rings have no damage or cracks	
			8. User manual User manual is available (User manuals can be downloaded from www.m	alift com)

### Approved by:

### Check all 8 points in checklist. If one or more of the points has a failure, the sling is not approved for use.



# **Periodic Inspection record**

Sling Model: (Ex. Tollet lowback)	Art.No:	Size:	Serial no.:	ок	
					t
				+	$^+$
				+	$^{+}$
					t
		M 2XX XS XS			
				+	4
				—	4
				—	4
				+	4
				1	
				—	4
				+	_
				+	_
				—	4
				—	4
				+	-
				+	1
				1	
				+	1
				1	
				+	1
				1	
				+	1
				1	
				1	1
				1	
				1	1
				1	
				1	1
				1	1

Equcation



Etac AS,

Destroy and discard damaged and old slings that does not pass inspection.



# Recording of periodic inspection on the sling

- Some manufacturer have a possibility to record the performed periodic inspection on the sling.
- Also a possibility to record when the sling is put into service for the first time







Considerations to make when choosing a sling

- The medical condition of the client
  - Does the client have muscle tonus?
  - Is the client able to sit without support?
  - Is the client able to sit at all?
  - Does the client have head control?
  - Does the client have involuntary movements?
  - Is the client able to understand and take command?
  - Does the patient have any decubitus?
  - Breathing affected?
  - Does the patient have leg amputees?
  - Pain issues?
  - How is the sensibility of the patient?



Education



Considerations to make when choosing a sling

- What kind of lift will be used? Lifting height?
  - Over-head/ceiling lift
  - Floor-based lift
  - Sit to stand lift
- What environment the sling has to be used
  - Mainly in a dry environment
  - For bathing
  - In a toilet situation
  - in an area with high demands for infection control
- Which transfers will be done?
  - From sitting to sitting
  - From lying to sitting
  - Lifting to and from the floor





Considerations to make when choosing a sling

- The size of the sling
  - Consequences when the sling is too small
    - The sling is squeezing the client and thereby causing pain
    - The leg rest are too short
    - The backrest is too short
    - The lifting straps are too short and the spreader bar is in the face of the client
  - Consequences when the sling is too big
    - There is an increased tendency for the client to slide through the sling
    - The sling gives less support
    - The lifting straps are too long and the lifting height is compromised



Considerations to make when choosing a sling

- The spreader bar 2, 3 or 4 point
  - Advantages with a 3 or 4 point spreader bar
    - Stability
    - More open space for the client
    - Highly recommended for bariatric users
    - Easy to reach the hooks when lifting from a lying position
- The spreader bar width
  - The width of the spreader bar can effect how much stability and support the sling gives the client
  - If the patient has pain over the shoulders a wider spreader bar can be helpful to avoid pressure.
  - During gait training, a spreader bar with a width of the shoulders gives a better support of the walking vest.



Education



Considerations to make when choosing a sling

• Tilting the patient in the sling

### In a reclined position

- ✓ Difficult to position the patient correct in a straight chair or toilet
- Tilt in space wheelchair is recommendable
- $\checkmark$  Less pressure on the legs More on the back
- ✓ For patient with severe reduced upper body stability
- Less tendency to fall forward out of the sling

Education

Recommendable to use a HighBack sling

### In an slightly reclined position

- Increased risk of not positioning the patient correct in a straight chair or toilet
- Less pressure on the legs More on the back
- For patient with reduced upper body stability
- ✓ The patient feels more secure
- Less tendency to fall forward out of the sling

### In an upright position

- ✓ Easier to access a straight chair or toilet
- $\checkmark$  A lot of pressure on the legs less on the back
- The patient must have good upper body stability
- ✓ The patient could feel more insecure
- Tendency to fall forward out of the sling



### References

BS EN ISO 10535:2006 Hoists for the transfer of disabled persons - Requirements and test methods  $2^{nd}$  edition 2006-12-15

Hop 6 The guide to the Handling of People a systems approach 6<sup>th</sup> edition Published by the National Back Pain Assosiation

Arbeta med personlyft – skrift, Reifeldt, Skebäck, Published by HjälpmedelsCenter Väst 2009





### Thank you for your attention!



### ole.thy@etac.com



